

UE Medical Airway Management Platform

Advanced Visualization | Clinical Confidence | Intelligent Airway Solutions

Empower Healthcare Professionals to Save Lives

AGENDA

Airway Management Challenges

Excellence Medical Care (TamiozMed) in Brief

UE Medical Company Overview

UE Integrated Airway Management Platform

UEScope Video Laryngoscope

SaCoVLM™ Video Laryngeal Mask

Flexible Bronchoscope: Single-use & Reusable

Why UE Medical?

Airway Management Challenges

Difficult Airway

Unexpected difficult anatomy can rapidly turn a routine intubation into a life-threatening emergency.

Time-Critical Procedures

Delays in securing the airway increase the risk of hypoxia, cardiac arrest, and poor patient outcomes.

First-Pass Success

Each additional intubation attempt significantly increases complications, trauma, and infection risk.

Limited Visualization

Traditional direct laryngoscopy often provides a restricted view of the glottis, especially in difficult airways.

Operator Variability


Clinical success depends heavily on the experience and skill of the clinician performing the procedure.

Clinical & Economic Impact

Failed or prolonged intubation leads to higher complication rates, longer hospital stays, and increased healthcare costs.

Airway Management Challenges

What if every clinician could achieve a clearer view, greater confidence, and a higher first-pass success rate?



Excellence Medical Care (TamiozMed) in Brief

A close-up photograph of a hand holding a blue paper airplane, symbolizing vision and forward movement. The hand is positioned on the left side of the frame, with the fingers gently gripping the wings of the plane. The background is a soft, out-of-focus blue sky with light clouds.

To inspire a healthier future by being the leading provider of innovative healthcare solutions across Saudi Arabia and the Gulf Region.

Vision

Excellence Medical Care (TamiozMed) in Brief

Driven by excellence and innovation, we partner with healthcare providers to deliver world-class medical technologies that create safer procedures, greater clinical confidence, and better patient care.

Mission

Excellence Medical Care (TamiozMed) in Brief

Values

- **Innovation:** Continuously advancing healthcare technologies.
- **Quality:** Delivering excellence without compromise.
- **Integrity:** Building trust through transparency and accountability.
- **Customer Success:** Creating long-term value for healthcare partners.
- **Patient First:** Supporting better clinical outcomes in every solution.

Excellence Medical Care (TamiozMed) in Brief

Ergomy
CARBON



Long life
SURGICAL INDUSTRIES



OMECA
Arab Medical Equipment
Company



Ortholand®

UE Medical Company Overview




About Us

Established in 2010, UE Medical has emerged as a leading provider of visualized medical devices, dedicated to serving clinical specialties including anesthesiology, pulmonology and ENT. Supported by a global network that includes 4 R&D centers, 3 sales centers, and 2 production bases, equipped with advanced EMC Lab, Optical Lab, EB Radiation Sterilization Facility and Automatic Machine, UE medical aims to play a significant role in the world's medical industry.

UE Medical ranks No.1 for airway management products in China, and top 3 providers of single-use video laryngoscopes in EMS across the USA, our products are trusted in over 11,000 hospitals, including hospitals in the USA, Germany, Spain, South Korea, and Japan, etc.

UE Medical Company Overview



UE Medical Company Overview

Founded in
2010
Xianju, Zhejiang

Over
580+
employees all over the
world

Cooperates with
80%
of Top tier
hospitals in China



Serving
70+
countries and regions

150+
R&D engineers
worldwide

11,000+
Hospitals
world-wide


TOP 3
Single-use
video laryngoscope
in Emergency in the
U.S.

4
R&D centers

2
Service centers


3
Sales centers

UE Medical Company Overview




UE Integrated Airway Management Platform

Scene Sharing for Enhanced Instruction



Data Sharing for Mobile Teaching



UE Integrated Airway Management Platform




One Monitor Compatible with All UE Products
Meet the comprehensive visualization needs of airway management.

Data management, process traceability

64G storage, supporting the collection, storage, query and sharing of anesthesia airway management data
Helping to trace the entire clinical diagnosis and treatment process



UE Integrated Airway Management Platform



Multi-Product Compatibility & Split-Screen Display


Supports simultaneous use of multiple products with a split-screen display, providing more comprehensive access to information.

Integrated Trolley for Efficiency & Convenience

Supports the storage and management of common airway management devices, consumables, and accessories.



UE Integrated Airway Management Platform




Versatile Applications for Various Scenarios

Supports in-hospital and out-of-hospital anesthesia airway consultations. Ideal for teaching, training, and surgical academic observations.

UE Integrated Airway Management Platform



UE Integrated Airway Management Platform



Video Laryngoscope (UED-D)



Video Laryngoscope (UED-A)



Video Laryngeal Mask




Video Laryngoscope (UED-C)



Disposable Video Bronchoscope


UEScope Video Laryngoscope



SaCoVLM™ Video Laryngeal Mask



SaCoVLM™ Video Laryngeal Mask



Flexible Bronchoscope: Single-use & Reusable



Flexible Bronchoscope: Single-use & Reusable

Monitor




Features

Large Bending Angle

Good for airway inspection

Operation part

Lightweight, portable and comfortable holding
Various model for choice



Camera

Clear image

Grab and go

Independent packaging, ready to use
No risk of cross-infection
Various specifications
Use in different scenarios

Health economy

Low input cost
No need for reprocessing equipment

Emergency Use

Portable structure,
Start the device in 3 seconds

Model & Recommended Uses

Insertion tube with extra small outer diameter
No working channel



PL-F260: 2.6mm
TIC-SD-I: 2.8mm
TIC-I1: 3.8mm

Uses:

1. inspection
2. pediatric airway management
3. check the placement of DLT
4. check the placement of endobronchial blocker

Insertion tube with small outer diameter
Working channel for drug administration



PL-F320: 3.1mm/1.2mm
PL-F400: 3.9mm/1.5mm
TIC-SD-III: 4.5mm/1.5mm

Uses:

1. endotracheal drug administration & inspection
2. tracheal intubation
3. check the placement of DLT

Insertion tube with large outer diameter
Working channel for suctioning



PL-F520: 5.2mm/2.6mm
TIC-I3: 5.2mm/2.2mm
TIC-I4: 5.5mm/2.4mm

Uses:

1. suction sputum
2. BAL procedures
3. remove foreign objects
4. sample lung issue

Why UE Medical?

- **Proven Market Leadership:** Trusted by thousands of hospitals and clinicians worldwide.
- **Comprehensive Airway Platform:** An integrated portfolio for video laryngoscopy, intubation stylets, bronchoscopy, and laryngeal masks.
- **Clinical Confidence:** Designed to enhance visualization, improve first-pass success, and support patient safety.
- **Continuous Innovation:** Advanced technologies developed to meet evolving clinical needs.
- **Outstanding Value:** High-quality solutions that deliver long-term clinical and economic benefits.
- **Reliable Partnership:** Dedicated clinical education, technical support, and responsive after-sales service.

THANK YOU

